

SAFETY CANS & CONTAINERS

for Laboratories



Tilt cradle allows can to swing forward for dispensing while holding can securely on bench top.

Safety Dispensing Cans

With all the benefits of our classic Type I safety cans (see page 4), these self-closing and leakproof cans come with faucets to make them the choice of lab professionals for safe dispensing of hazardous liquids into small flasks and beakers. Fill spouts have self-closing pressure-relief caps and stainless steel flame arresters. Self-closing brass faucets provide trouble-free dispensing. Cans fitted with the 08540 faucet offer more controlled dispensing, while faucet 08902 provides a greater flow.

Polyethylene Safety Dispensing Can has superb corrosion protection, and a unique current-carrying carbon insert embedded into the container allows for proper grounding.

Steel & Poly Shelf Cans have a faucet mounted 2-in (51-mm) above the bottom to eliminate the need to lift the cans during dispensing.

Steel Tilt Cans have a faucet on top and the can rests in a sturdy metal cradle. Cradle base measures 13-13/16 x 14-in (351-mm x 356-mm).

All styles are FM approved.



Store & secure chemicals . . .

For security against unauthorized use, store safety cans in lockable Sure-Grip® EX safety cabinets – see pages 45-52.



7150250



14545Z



7150156

7150140



Description and Capacity	Can Dimensions OD x H	Approv/Lstg Regulation	Red Can Model No	Yellow Can Model No	Ship Wt lb/kg
2.5 gallon (9.5 L) steel safety shelf can with bottom faucet 08540	11.75 x 11.5 in 298 x 292 mm	FM	7225140	—	9/4
2.5 gallon (9.5 L) steel safety shelf can with bottom faucet 08902	11.75 x 16.875 in 298 x 429 mm		7225150Z	7225250	
5 gallon (19 L) steel safety shelf can with bottom faucet 08540			7150140	7150240	
5 gallon (19 L) steel safety shelf can with bottom faucet 08902			7150150Z	7150250	
5 gallon (19 L) tilt steel safety can with top faucet 08540 and stand	11.75 x 15.875 in 298 x 403 mm	FM, UL	7150146	—	13/6
5 gallon (19 L) tilt steel safety can with top faucet 08902 and stand	7150156		—		
5 gallon (19 L) poly safety shelf can with self-close brass faucet 08540	11.5 x 17 in 292 x 432 mm	FM	14545Z	—	11/5