

1. Part Number: 8085A-100S

2. Photograph:



3. Product Information:

Description:	Disposable Sack Sampler – 400mm - Sterile
Use:	Single use, disposable powder sampler
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004
Method of Construction:	Single piece construction
Overall Length:	400mm
Diameter of Outer Tube:	25mm
Nominal Weight of Sampler:	53g
Quantity per Box:	100
Individually Bagged?	Yes (heat sealed PE bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	2 years from the date of irradiation

4. Change Control:

Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill, Warwickshire, B46 1JH, UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21