



Mini-100 Orbital-Genie™

Compact slow orbital shaking on a non-slip platform

Non-slip, platform designed with spill-proof edges is ideal for gentle shaking applications such as Western Blots, with no additional accessories necessary. Variable shaking speed control up to 100RPM.



	100 RPM with Ultra High Torque
	Super Compact Design
	Variable Speed
	Gentle Orbital Shaking

Accessories

A non-slip tray that can be stacked onto the Mini-100 Orbital-Genie's existing platform tray, increasing the capacity for mixing by stacking up two trays. No tools required, stays in place, simply lift-to-remove.



Add stacking trays

Double or triple the capacity of your Mini-100



Specifications

	Mini-100 Orbital-Genie
Speed (RPM)	20 – 100 RPM
Orbit	19mm
Weight	3.5 Kg (7.8 lbs.)
Dimensions	350 x 350mm (14 x 14in)
Platform	290 x 290mm (11.4 x 11.4in)
Platform Weight Capacity	3 Kg (6.6 lbs)
Markings	CE

Mini Orbital-Genie Ordering Information

	Catalog No.
Voltage & Plug configuration	Mini-100 Orbital-Genie
120V - U.S Standard	SI-M100
230V - No Plug	SI-M101
230V - European Plug	SI-M102
230V - British Plug	SI-M103
230V - Swiss Plug	SI-M104
240V - Aust. Plug	SI-M106
100V	SI-M105

Accessory Ordering Information

Catalog No.	Description
SI-M150	Stacking Tray, for the Mini-100 Orbital-Genie



80 Orville Drive,
Suite 102,
Bohemia,
New York 11716 USA



Phone: (631) 567-4700
Toll Free: (888) 850-6208



Customer Service Email:
custsvc@scientificindustries.com



Fax: (631) 567-5896

