

# VORNADO™

## Miniature Vortex Mixer



- + Powerful vortexing for tubes up to 50ml
- + Insta-Touch™ pressure activated operation
- + "Mini" sized, less than 4 in. wide



Extremely Compact



Insta-Touch™ Pressure Activation

**Powerful** is the best word to describe this **mini** vortexer. With a 4mm orbit and fixed speed of 2800rpm, the Vornado instantly vortexes even the largest samples, including nearly full 50 ml tubes. Despite its powerful motor, the Vornado has a footprint of less than 4x4 inches, allowing it to fit on even the most crowded bench.

Insta-Touch operation means that the motor is automatically activated with slight pressure on the center cup head. Operation is smooth, quiet and efficient. When pressure is removed, operation immediately ceases.

The Vornado's unique head design prevents liquids from entering the housing, prolonging motor life. The unit is available with five different cup head colors.

### Technical Data:

Speed:	2800 rpm
Motor Activation:	Touch
Orbit:	4mm
Operating Temp. Range:	+4°C to +65°C
Dimensions:	9.4 x 9.9 x 6.6 cm
	3.7 x 3.9 x 2.6 in.
Weight:	0.4 kg / 1 lb
Electrical:	100-240V, 50-60 Hz, 0.7A
Warranty:	2 years



BV101

BV101-B



BV101-G

BV101-P

BV101-R

### Ordering Information:

BV101*	Vornado™ Miniature Vortex Mixer, grey cup head
BV101-B*	Vornado™ Miniature Vortex Mixer, blue cup head
BV101-G*	Vornado™ Miniature Vortex Mixer, green cup head
BV101-P*	Vornado™ Miniature Vortex Mixer, purple cup head
BV101-R*	Vornado™ Miniature Vortex Mixer, red cup head

\* Includes US Plug. For EU plug, please add (-E)

**Benchmark**  
Scientific

**Benchmark Scientific Inc.**  
PO Box 709 Edison, NJ 08818  
Phone: 908-769-5555  
Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)